


New
Puts an end to Mortellaro
and other hoof diseases



Patent pending

Vink
Livestock equipment

HYDROMATIC FOOT BATH from VINK


This walk-through foot bath consists of two adjacent channels through which your medication or disinfectant is pumped at great speed. When the cows walk through the channels, the liquid is sprayed forcefully against the feet and heel removing contamination from wounds and the interdigital space. The powerful jets spray the medication into the hard-to-reach areas between the cow's toes and behind her heel.




The tray for treating cow's feet uses two divided channels which helps the cows to walk slowly through the 3 metre long bath. This process allows the treatment sufficient time to penetrate her feet thoroughly.

There is a small holding tank connected to the channels which contains 15 litres (4 gallons) of liquid. This holding tank is constantly replenished by a reservoir located elsewhere in the barn (see diagram opposite). The unit is activated when the cow enters the foot bath. As the cow walks through the foot bath and receives her treatment, the water is filtered with a sieve and brush to separate the solid contaminants from the liquid and recycles the water back into the small holding tank which is continuously topped up to 15 litres (4 gallons) from the reservoir, which is already filled with water and a disinfectant or medication. This maintains the quality of the treatment from the first cow until the last cow.

The cows quickly become familiar with the tray, however it is recommended that cows walk through the dry foot bath twice before the spray mechanism is put into operation.





The system contains about 15 litre of liquid, which is pumped around as soon as the sensor detects the cow. Approx. 3/4 litre per cow spilled liquid will be automatically refilled from the reservoir

The advantages:

- Powerful jets spray the liquid against the cow's foot and heel for deep cleaning in the critical interdigital space
- The design and length of the channels help the cows to walk slowly, therefore achieving intensive cleaning and penetration
- The filtration techniques along with continuously adding new liquid to the holding tank preserves the quality of the treatment from the first cow to the last cow.
- **No manual work** - the system cleans and fills itself automatically
- The system is solidly constructed entirely of stainless steel
- All substances that are soluble in water can be used
- Water consumption is low – Approx. 4 litre (1 gallon) per 5 cows.



The two divided channels help the cows to walk slowly and carefully through the 3 metre long bath. This process allows the treatment sufficient time to penetrate her feet thoroughly.



When the cows walk through the channels, the liquid is sprayed forcefully against the feet and heel removing contamination from wounds and the interdigital space.

Position:

Permanently installed after exit from the milking parlour.

The tray for treating cow's feet must be installed in such a place that when it is not in use, the cows can exit around it

Use:

For curative treatment:

3 to 7 x per week

Preventative, treatment:

1x per week

After 6 weeks of curative treatment, cows will be free of infectious foot disorders.

Water consumption:

Approx. 4 litre (1 gallon) per 5 cows.

Remember, the tray for treating cow's feet is automatically cleaned and the water recycled and refilled continuously.

Material:

The 3.2 metre foot bath is completely stainless steel

The tray for treating cow's feet by Vink is an ideal method for treating and preventing infectious hoof disorders.

Vink Livestock equipment

The logo for Vink, featuring the word "Vink" in a bold, blue, sans-serif font. The letter "i" has a red dot above it.